1A400

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
Assembling Instructions

1- Description and functions

Dyna C foot includes three main components. Two composite parts: a C-form malleolus that allow energy storage and return, and a sole providing flexibility on the side, which ensure adaptation to the ground. A polyurethane cosmetic shell which transfers the mechanical stresses between the ground or the shoe and the composite structure. Thanks to these components, the foot behaviour is flexible and dynamic.

For patients who are active, or very active (children for example), or reach the weight and activity limits of the module selection guide, a stop can be used to limit warping of the C-malleolus. This stop comes with the foot and must be installed by the prothesist.

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Dyna C foot is designed exclusively for lower limb prostheses. It is appropriate for active patients up to 125kg, or very active patients up to 105kg.


Dyna C foot is designed for a 10mm heel. It includes a male pyramid that allows for the use of standard connectors available on the market. Thanks to these connectors, the heel height can be easily modified from 0 to 20mm.

Module selection guide


Activité	Patient weight (kg)				
	50 to 60	60 to 75	75 to 90	90 to 105	105 to 125
Moderate	1	1	2	3	4
Normal	1	2	3	4	5
High	2	3	4	5	N/A

3- Montage (sans le patient)

A) Gluing of the stop (if necessary): Remove the foot from its cosmetic shell, then take the sock off. Degrease the faces that will be glued, let them dry, then apply epoxy glue on both faces.

Put the stop in contact with the nut (1), then position it by pushing it upwards (2) and toward the nut (3). Let the glue dry, then put the sock and the cosmetic shell again.

B) Assembling: Place under the foot a wedge of the same height as the shoe heel. Fix the female pyramid connector to the foot taking care that its axis be perpendicular to the ground (fig 1).


Fig. 1

4- Alignment (with the patient):

Place the foot so that the load line goes through the posterior third of the foot (fig. 2).



Fig. 2

Notes:

- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if the heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems too hard, or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

5- Finishing:

To avoid rubbing noises, pull the sock up to the ankle connector and fix it with a fine sting for example. The malleolar clips that ensure connection to the cosmetic foam allow to dismantle the foot easily and rapidly.

	Foot size (cm)	22 to 25	26 to 27	28 to 29
	Product Nr	1A40013	1A40014	1A40015


Installation: Cut the foam to the desired length, position the foam on the foot and draw the malleolus contour. Inside the foam mill a space that will receive the selected connection part and the foot part protruding from the shell, and at the same time take care that thickness be sufficient to allow for gluing. Position the malleolar clip carefully and clip it into cosmesis. Clean contact surfaces with isopropyl alcohol and glue the malleolar clip on the foam with XC090 néoprene adhesive. After drying, finish the outer shape as usualy.

6- Maintenance:

The cosmetic shell, the protective sock and the malleolar clip are the only foot components that may be worn out and need to be replaced. The cosmesis can be removed and plug again with 1A40019 insert horn.

Spare cosmesis 1A40020-D22 including the cosmetic shell, the sock and the malleolar clip.

Side  Size in cm

7- Recycling:

The different components of the foot are special wastes and must be treated accordingly.

Part for the patient

ADVICE FOR USE, MAINTENANCE AND SAFETY

Clean the cosmetic cover with a slightly damp cloth or sponge.

Use temperature : - 20°C / +60 °C.

Do not let foreign bodies or dust inside the foot.

If the foot has been dipped or sprinkled, empty the remaining liquid and dry the foot rapidly.

A use of the foot that does not adhere to your prosthetist's recommendations can damage foot components.

In case of an abnormal behaviour of your foot, or changes in your foot features, or if the foot sustained a hard shock, contact your prosthetist.

PROTEOR
6 rue de la Redoute - 21850 Saint-Apollinaire - France
Tél. 33 (0) 3 80 78 42 10 - Fax 33 (0) 3 80 78 42 15
ht.orthopedie@proteor.com - www.proteor.com

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